


Platinum Labware



Precious metals and catalysts

Browse products >



Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Саранск (8342)22-96-24
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(7172)727-132

Киргизия +996(312)96-26-47

<https://aesar.nt-rt.ru/> || arj@nt-rt.ru

Care and Use of your Platinum Labware



Platinum ware can be heated in many atmospheres commonly employed for laboratory experimentation such as:

Air, nitrogen and oxygen
Oxides of nitrogen
Bromine and iodine vapors
Hydrogen
Carbon dioxide

Platinum ware is suitable for numerous chemical research applications such as ignitions, fusions, evaporations /hand glass production. These are usually accomplished in oxidizing atmospheres.

37995

37995 Cleaning Kit for Laboratory Crucibles

Platinum Crucible for Leco FX-100 & FX-200 Automatic Fluxer with Supporting Ring



Platinum Labware

[Request a Quote >](#)



46790

Platinum crucible for Leco FX-100 & FX-200 auto fluxer w/sup ring, Top Dia 35mm, Bot Dia 21mm, Ht 39mm, Bt 0.60mm



46702

Platinum crucible for Leco FX-100 & FX-200 auto fluxer w/sup ring, Top Dia 37mm, Bot Dia 20mm, Ht 35mm, Bt 0.60mm

Platinum Lids for Micro Crucibles



Platinum Labware

[Request a Quote >](#)



46686 Platinum lid for micro crucible, Dia 12mm, fits 46413



46993 Platinum lid for micro crucible, Dia 17mm, fits 46762



46349 Platinum lid for micro crucible, Dia 20.5mm, fits 46472



46776 Platinum lid for micro crucible, Dia 26mm, fits 46785

Platinum Standard Form Crucibles with Reinforced Rim



High purity metals and materials
Fundamentals for research



[Discover our range >](#)

Platinum Labware



[Request a Quote >](#)

46956 Platinum crucible with reinforced rim, Top Dia 26mm, Bot Dia 15.5mm, Ht 29mm, Base Thickness 0.25mm, Capacity 10mL

46591 Platinum crucible with reinforced rim, Top Dia 29mm, Bot Dia 17.5mm, Ht 24mm, Base Thickness 0.25mm, Capacity 10mL

46468	Platinum crucible with reinforced rim, Top Dia 29mm, Bot Dia 17.5mm, Ht 33mm, Base Thickness 0.25mm, Capacity 15mL
46577	Platinum crucible with reinforced rim, Top Dia 32.5mm, Bot Dia 18mm, Ht 36mm, Base Thickness 0.3mm, Capacity 20mL
46489	Platinum crucible with reinforced rim, Top Dia 32mm, Bot Dia 19.5mm, Ht 27mm, Base Thickness 0.25mm, Capacity 15mL
46353	Platinum crucible with reinforced rim, Top Dia 35mm, Bot Dia 21mm, Ht 39mm, Base Thickness 0.3mm, Capacity 25mL
46580	Platinum crucible with reinforced rim, Top Dia 37.5mm, Bot Dia 22.5mm, Ht 41mm, Base Thickness 0.3mm, Capacity 30mL
46894	Platinum crucible with reinforced rim, Top Dia 37mm, Bot Dia 22.5mm, Ht 31mm, Base Thickness 0.25mm, Capacity 20mL
46965	Platinum crucible with reinforced rim, Top Dia 39mm, Bot Dia 23.5mm, Ht 43mm, Base Thickness 0.3mm, Capacity 35mL
46958	Platinum crucible with reinforced rim, Top Dia 41.5mm, Bot Dia 25mm, Ht 46mm, Base Thickness 0.3mm, Capacity 40mL
46618	Platinum crucible with reinforced rim, Top Dia 41mm, Bot Dia 25mm, Ht 34mm, Base Thickness 0.3mm, Capacity 25mL
46722	Platinum crucible with reinforced rim, Top Dia 43.5mm, Bot Dia 26.5mm, Ht 36mm, Base Thickness 0.3mm, Capacity 35mL
46561	Platinum crucible with reinforced rim, Top Dia 44mm, Bot Dia 26.5mm, Ht 49mm, Base Thickness 0.38mm, Capacity 50mL
46719	Platinum crucible with reinforced rim, Top Dia 49mm, Bot Dia 30mm, Ht 41mm, Base Thickness 0.36mm, Capacity 50mL
46517	Platinum crucible with reinforced rim, Top Dia 50mm, Bot Dia 30mm, Ht 55mm, Base Thickness 0.41mm, Capacity 75mL
46633	Platinum crucible with reinforced rim, Top Dia 55mm, Bot Dia 33.5mm, Ht 46mm, Base Thickness 0.4mm, Capacity 70mL
46883	Platinum crucible with reinforced rim, Top Dia 56mm, Bot Dia 33mm, Ht 62mm, Base Thickness 0.48mm, Capacity 100mL
46737	Platinum crucible with reinforced rim, Top Dia 59mm, Bot Dia 35mm, Ht 64mm, Base Thickness 0.48mm, Capacity 125mL
46336	Platinum crucible with reinforced rim, Top Dia 62mm, Bot Dia 37mm, Ht 68mm, Base Thickness 0.48mm, Capacity 150mL
46635	Platinum crucible with reinforced rim, Top Dia 65mm, Bot Dia 39mm, Ht 71mm, Base Thickness 0.48mm, Capacity 175mL

Platinum Micro Crucibles



Platinum Labware

[Request a Quote >](#)



46413 Platinum micro crucible, Top OD 12mm, Bot Dia 7.2mm, Ht 13mm, Capacity 1mL



46762 Platinum micro crucible, Top OD 17mm, Bot Dia 10mm, Ht 19mm, Capacity 3mL



46472 Platinum micro crucible, Top OD 20.5mm, Bot Dia 12mm, Ht 22mm, Capacity 5mL



46785 Platinum micro crucible, Top OD 26mm, Bot Dia 15.5mm, Ht 28mm, Capacity 10mL

Platinum Standard Lids for Form Crucibles



High purity metals and materials
Fundamentals for research



[Discover our range >](#)

Platinum Labware



[Request a Quote >](#)



46961

Platinum Lid for Standard form crucible, Dia 19.5mm, fits Stock # 46916



46535

Platinum Lid for Standard form crucible, Dia 22.5mm, fits Stock # 46377

46337	Platinum Lid for Standard form crucible, Dia 25.5mm, fits Stock # 46800
46979	Platinum Lid for Standard form crucible, Dia 26mm, fits Stock #s 46800 & 46956
46678	Platinum Lid for standard form crucible, Dia 29.5mm, fits Stock # 46484
46401	Platinum Lid for Standard form crucible, Dia 29mm, fits Stock #s 46825 & 46468
46391	Platinum Lid for Standard form crucible, Dia 29mm, fits Stock #s 46896 & 46591
46739	Platinum Lid for Standard form crucible, Dia 32.5mm, fits Stock #s 46520 & 46577
46407	Platinum Lid for Standard form crucible, Dia 32mm, fits Stock #s 46371 & 46489
46395	Platinum Lid for Standard form crucible, Dia 34.5mm, fits Stock #s 46375 & 46345
46364	Platinum Lid for Standard form crucible, Dia 35mm, fits Stock #s 46870 & 46353
46732	Platinum Lid for Standard form crucible, Dia 37.5mm, fits Stock #s 46813 & 46580
46514	Platinum Lid for Standard form crucible, Dia 37mm, fits Stock #s 46498 & 46894
46402	Platinum Lid for Standard form crucible, Dia 39mm, fits Stock #s 46649 & 46965
46355	Platinum Lid for Standard form crucible, Dia 41.5mm, fits Stock #s 46911 & 46958
46626	Platinum Lid for Standard form crucible, Dia 41mm, fits Stock #s 46578 & 46618
46710	Platinum Lid for Standard form crucible, Dia 43.5mm, fits Stock #s 46477 & 46722
46913	Platinum Lid for Standard form crucible, Dia 44.5mm, fits Stock # 46436
46868	Platinum Lid for Standard form crucible, Dia 44mm, fits Stock #s 46931 & 46561
46367	Platinum Lid for Standard form crucible, Dia 49.5mm, fits Stock # 46607
46554	Platinum Lid for Standard form crucible, Dia 49mm, fits Stock #s 46717 & 46719

46744	Platinum Lid for Standard form crucible, Dia 50mm, fits Stock #s 46420 & 46517
46655	Platinum Lid for Standard form crucible, Dia 54.5mm, fits Stock # 46670
46725	Platinum Lid for Standard form crucible, Dia 55mm, fits Stock #s 46793 & 46633
46567	Platinum Lid for Standard form crucible, Dia 56mm, fits Stock #s 46509 & 46737
46427	Platinum Lid for Standard form crucible, Dia 59.5mm, fits Stock # 46521
46843	Platinum Lid for Standard form crucible, Dia 59mm, fits Stock #s 46857 & 46737
46382	Platinum Lid for Standard form crucible, Dia 62mm, fits Stock #s 46631 & 46336
46735	Platinum Lid for Standard form crucible, Dia 65mm, fits Stock #s 46500 & 46635
46694	Platinum Lid for Standard form crucible, Dia 68mm, fits Stock #s 46465 & 46656
46552	Platinum Lid for Standard form crucible, fits Stock # 46964

Platinum Roasting Crucible with Lid



High purity metals and materials
Fundamentals for research



[Discover our range >](#)

Platinum Labware



[Request a Quote >](#)



46768

Platinum roasting crucible with lid, Base Thickness 0.25mm, Capacity 30mL

Platinum Straight Wall Crucible Lids



High purity metals and materials
Fundamentals for research



[Discover our range >](#)

Platinum Labware



[Request a Quote >](#)



47313 Platinum Lid for straight wall crucible, fits Stock # 46432

Platinum Standard Form Crucibles



High purity metals and materials
Fundamentals for research



[Discover our range >](#)

Platinum Labware



[Request a Quote >](#)








46916













Platinum standard form crucible, Top OD 20.5mm, Bot Dia 12mm, Ht 22mm, 5mL



46377

Platinum standard form crucible, Top OD 23mm, Bot Dia 14mm, Ht 27mm, Capacity 7mL

	46751	Platinum standard form crucible, Top OD 26mm, Bot Dia 15.5mm, Ht 29mm, Base Thickness 0.25mm, Capacity 10mL
	46800	Platinum standard form crucible, Top OD 26mm, Bot Dia 17mm, Ht 30mm, Capacity 10mL
	46974	Platinum standard form crucible, Top OD 28mm, Bot Dia 18mm, Ht 32mm, Capacity 15mL
	46896	Platinum standard form crucible, Top OD 29mm, Bot Dia 17.5mm, Ht 24mm, Base Thickness 0.25mm, Capacity 10mL
	46825	Platinum standard form crucible, Top OD 29mm, Bot Dia 17.5mm, Ht 33mm, Base Thickness 0.25mm, Capacity 15mL
	46484	Platinum standard form crucible, Top OD 30mm, Bot Dia 19mm, Ht 35mm, Capacity 20mL
	46520	Platinum standard form crucible, Top OD 32.5mm, Bot Dia 18mm, Ht 36mm, Base Thickness 0.3mm, Capacity 20mL
	46371	Platinum standard form crucible, Top OD 32mm, Bot Dia 19.5mm, Ht 27mm, Base Thickness 0.25mm, Capacity 15mL
	46870	Platinum standard form crucible, Top OD 35mm, Bot Dia 21mm, Ht 39mm, Base Thickness 0.3mm, Capacity 25mL
	46375	Platinum standard form crucible, Top OD 35mm, Bot Dia 23mm, Ht 37mm, Capacity 25mL
	46343	Platinum standard form crucible, Top OD 35mm, Bot Dia 23mm, Ht 40mm, Capacity 30mL
	46813	Platinum standard form crucible, Top OD 37.5mm, Bot Dia 22.5mm, Ht 41mm, Base Thickness 0.3mm, Capacity 30mL
	46498	Platinum standard form crucible, Top OD 37mm, Bot Dia 22.5mm, Ht 31mm, Base Thickness 0.25mm, Capacity 20mL
	46649	Platinum standard form crucible, Top OD 39mm, Bot Dia 23.5mm, Ht 43mm, Base Thickness 0.3mm, Capacity 35mL
	46911	Platinum standard form crucible, Top OD 41.5mm, Bot Dia 25mm, Ht 46mm, Base Thickness 0.3mm, Capacity 40mL
	46964	Platinum standard form crucible, Top OD 41.5mm, Bot Dia 25mm, Ht 46mm, Capacity 40 ml
	46578	Platinum standard form crucible, Top OD 41mm, Bot Dia 25mm, Ht 34mm, Base Thickness 0.3mm, Capacity 25mL
	46477	Platinum standard form crucible, Top OD 43.5mm, Bot Dia 26.5mm, Ht 36mm, Base Thickness 0.3mm, Capacity 35mL
	46931	Platinum standard form crucible, Top OD 44mm, Bot Dia 26.5mm, Ht 49mm, Base Thickness 0.38mm, Capacity 50mL

	46436	Platinum standard form crucible, Top OD 45mm, Bot Dia 28mm, Ht 50mm, Capacity 50mL
	46717	Platinum standard form crucible, Top OD 49mm, Bot Dia 30mm, Ht 41mm, Base Thickness 0.36mm, Capacity 50mL
	46420	Platinum standard form crucible, Top OD 50mm, Bot Dia 30mm, Ht 55mm, Base Thickness 0.41mm, Capacity 75mL
	46607	Platinum standard form crucible, Top OD 50mm, Bot Dia 32mm, Ht 55mm, Capacity 75mL
	46793	Platinum standard form crucible, Top OD 55mm, Bot Dia 33.5mm, Ht 46mm, Base Thickness 0.4mm, Capacity 70mL
	46670	Platinum standard form crucible, Top OD 55mm, Bot Dia 36mm, Ht 60mm, Capacity 100mL
	46509	Platinum standard form crucible, Top OD 56mm, Bot Dia 33mm, Ht 62mm, Base Thickness 0.48mm, Capacity 100mL
	46857	Platinum standard form crucible, Top OD 59mm, Bot Dia 35mm, Ht 64mm, Base Thickness 0.48mm, Capacity 125mL
	46521	Platinum standard form crucible, Top OD 60mm, Bot Dia 40mm, Ht 65mm, Capacity 125mL
	46631	Platinum standard form crucible, Top OD 62mm, Bot Dia 37mm, Ht 68mm, Base Thickness 0.48mm, Capacity 150mL
	46500	Platinum standard form crucible, Top OD 65mm, Bot Dia 39mm, Ht 71mm, Base Thickness 0.48mm, Capacity 175mL
	46465	Platinum standard form crucible, Top OD 68mm, Bot Dia 41mm, Ht 74mm, Base Thickness 0.48mm, Capacity 200mL

Platinum Straight Wall Crucibles



High purity metals and materials
Fundamentals for research



[Discover our range >](#)

Platinum Labware



[Request a Quote >](#)

46841

Platinum straight wall crucible, Top OD 18mm, Bot Dia 10mm, Ht 15mm, Base Thickness 0.2mm, Capacity 2mL

46771

Platinum straight wall crucible, Top OD 18mm, Bot Dia 10mm, Ht 15mm, Base Thickness 0.5mm, Capacity 2mL

	46411	Platinum straight wall crucible, Top OD 20mm, Bot Dia 15mm, Ht 20mm, Base Thickness 0.2mm, Capacity 4.5mL
	46705	Platinum straight wall crucible, Top OD 20mm, Bot Dia 15mm, Ht 20mm, Base Thickness 0.5mm, Capacity 4.5mL
	46827	Platinum straight wall crucible, Top OD 25mm, Bot Dia 17mm, Ht 30mm, Base Thickness 0.2mm, Capacity 10mL
	46480	Platinum straight wall crucible, Top OD 25mm, Bot Dia 17mm, Ht 30mm, Base Thickness 0.5mm, Capacity 10mL
	46808	Platinum straight wall crucible, Top OD 30mm, Bot Dia 20mm, Ht 55mm, Base Thickness 0.2mm, Capacity 25mL
	46432	Platinum straight wall crucible, Top OD 35mm, Bot Dia 21mm, Ht 39mm, Base Thickness 0.3mm, Capacity 25mL
	46572	Platinum straight wall crucible, Top OD 40mm, Bot Dia 30mm, Ht 50mm, Base Thickness 0.2mm, Capacity 45mL
	46920	Platinum straight wall crucible, Top OD 40mm, Bot Dia 30mm, Ht 50mm, Base Thickness 0.5mm, Capacity 45mL
	46505	Platinum straight wall crucible, Top OD 50mm, Bot Dia 40mm, Ht 50mm, Base Thickness 0.2mm, Capacity 75mL
	46887	Platinum straight wall crucible, Top OD 50mm, Bot Dia 40mm, Ht 50mm, Base Thickness 0.5mm, Capacity 75mL
	46743	Platinum straight wall crucible, Top OD 60mm, Bot Dia 50mm, Ht 55mm, Base Thickness 0.2mm, Capacity 130mL
	46992	Platinum straight wall crucible, Top OD 60mm, Bot Dia 50mm, Ht 55mm, Base Thickness 0.5mm, Capacity 130mL
	46698	Platinum straight wall crucible, Top OD 60mm, Bot Dia 55mm, Ht 35mm, Base Thickness 0.2mm, Capacity 90mL
	46829	Platinum straight wall crucible, Top OD 60mm, Bot Dia 55mm, Ht 35mm, Base Thickness 0.5mm, Capacity 90mL
	46369	Platinum straight wall crucible, Top OD 65mm, Bot Dia 55mm, Ht 75mm, Base Thickness 0.2mm, Capacity 210mL
	46755	Platinum straight wall crucible, Top OD 65mm, Bot Dia 55mm, Ht 75mm, Base Thickness 0.5mm, Capacity 210mL

Non-Wetting Platinum-Gold Apparatus



Flux fusions are widely employed on a routine basis for the preparation of samples for X-ray fluorescence and other instrumental and classical analysis techniques.

/n

Several types of automatic equipment are available for this work and specially designed crucibles and molds /nare produced for use with this equipment. Crucibles and molds for /nautomatic equipment such as Claisse, Perl X2, LECO, Herzog and other forms are also available. Normally, /nnon-wetting Pt 5% Au alloy is used for this work, although any alloy or material may be specified.

Non-Wetting Platinum 5% Gold Crucibles With Reinforced Rim



High purity metals and materials
Fundamentals for research



[Discover our range >](#)

Platinum Labware



[Request a Quote >](#)

46490

Non-wetting Pt 5% Au crucible with RF rim, Top Dia 26mm, BotDia 15.5mm, Ht 29mm, Base Thickness 0.25mm, Capacity 10mL

46988

Non-wetting Pt 5% Au crucible with RF rim, Top Dia 29mm, BotDia 17.5mm, Ht 24mm, Base Thickness 0.25mm, Capacity 10mL

46657	Non-wetting Pt 5% Au crucible with RF rim, Top Dia 29mm, BotDia 17.5mm, Ht 33mm, Base Thickness 0.25mm, Capacity 15mL
46558	Non-wetting Pt 5% Au crucible with RF rim, Top Dia 32.5mm, Bot Dia 18mm, Ht 36mm, Base Thickness 0.25mm, Capacity 20mL
46930	Non-wetting Pt 5% Au crucible with RF rim, Top Dia 32mm, BotDia 19.5mm, Ht 27mm, Base Thickness 0.30mm, Capacity 15mL
46741	Non-wetting Pt 5% Au crucible with RF rim, Top Dia 35mm, BotDia 21mm, Ht 39mm, Base Thickness 0.30mm, Capacity 25mL
46565	Non-wetting Pt 5% Au crucible with RF rim, Top Dia 37.5mm, Bot Dia 22.5mm, Ht 41mm, Base Thickness 0.30mm, Cap. 30mL
46659	Non-wetting Pt 5% Au crucible with RF rim, Top Dia 37mm, BotDia 22.5mm, Ht 31mm, Base Thickness 0.30mm, Capacity 20mL
46688	Non-wetting Pt 5% Au crucible with RF rim, Top Dia 39mm, BotDia 23.5mm, Ht 43mm, Base Thickness 0.30mm, Capacity 35mL
46963	Non-wetting Pt 5% Au crucible with RF rim, Top Dia 41.5mm, Bot Dia 25mm, Ht 46mm, Base Thickness 0.30mm, Capacity 40mL
46675	Non-wetting Pt 5% Au crucible with RF rim, Top Dia 41mm, BotDia 25mm, Ht 34mm, Base Thickness 0.30mm, Capacity 25mL
46495	Non-wetting Pt 5% Au crucible with RF rim, Top Dia 43.5mm, Bot Dia 26.5mm, Ht 36mm, Base Thickness 0.30mm, Cap. 35mL
46389	Non-wetting Pt 5% Au crucible with RF rim, Top Dia 44mm, BotDia 26.5mm, Ht 49mm, Base Thickness 0.38mm, Capacity 50mL
46982	Non-wetting Pt 5% Au crucible with RF rim, Top Dia 49mm, BotDia 30mm, Ht 41mm, Base Thickness 0.36mm, Capacity 50mL
46779	Non-wetting Pt 5% Au crucible with RF rim, Top Dia 50mm, BotDia 30mm, Ht 55mm, Base Thickness 0.41mm, Capacity 75mL
46459	Non-wetting Pt 5% Au crucible with RF rim, Top Dia 55mm, BotDia 33.5mm, Ht 46mm, Base Thickness 0.40mm, Capacity 70mL
46458	Non-wetting Pt 5% Au crucible with RF rim, Top Dia 56mm, BotDia 33mm, Ht 62mm, Base Thickness 0.48mm, Capacity 100mL
46804	Non-wetting Pt 5% Au crucible with RF rim, Top Dia 58mm, BotDia 35mm, Ht 64mm, Base Thickness 0.48mm, Capacity 125mL
46429	Non-wetting Pt 5% Au crucible with RF rim, Top Dia 62mm, BotDia 37mm, Ht 68mm, Base Thickness 0.48mm, Capacity 150mL
46850	Non-wetting Pt 5% Au crucible with RF rim, Top Dia 65mm, BotDia 39mm, Ht 71mm, Base Thickness 0.48mm, Capacity 175mL

Platinum Lids for Non-Wetting Crucibles



High purity metals and materials
Fundamentals for research



[Discover our range >](#)

Platinum Labware



[Request a Quote >](#)

46773

46773

Platinum Lid for non-wetting crucible, Dia 26mm, fits Stock #s 46376 & 46490

46620

46620

Platinum Lid for non-wetting crucible, Dia 29mm, fits Stock #s 46483 & 46988

46660	Platinum Lid for non-wetting crucible, Dia 29mm, fits Stock #s 46845 & 46657
46677	Platinum Lid for non-wetting crucible, Dia 32.5mm, fits Stock #s 46900 & 46558
46875	Platinum Lid for non-wetting crucible, Dia 32mm, fits Stock #s 46769 & 46930
46723	Platinum Lid for non-wetting crucible, Dia 35mm, fits Stock #s 46859 & 46741
46615	Platinum Lid for non-wetting crucible, Dia 37.5mm, fits Stock #s 46865 & 46565
46417	Platinum Lid for non-wetting crucible, Dia 37mm, fits Stock #s 46354 & 46659
46527	Platinum Lid for non-wetting crucible, Dia 39mm, fits Stock #s 46822 & 46688
46566	Platinum Lid for non-wetting crucible, Dia 41.5mm, fits Stock #s 46435 & 46963
46450	Platinum Lid for non-wetting crucible, Dia 41mm, fits Stock #s 46542 & 46675
46634	Platinum Lid for non-wetting crucible, Dia 43.5mm, fits Stock #s 46616 & 46495
46442	Platinum Lid for non-wetting crucible, Dia 44mm, fits Stock #s 46701 & 46389
46936	Platinum Lid for non-wetting crucible, Dia 49mm, fits Stock #s 46400 & 46982
46646	Platinum Lid for non-wetting crucible, Dia 50mm, fits Stock #s 46901 & 46779
46531	Platinum Lid for non-wetting crucible, Dia 55mm, fits Stock #s 46840 & 46459
46586	Platinum Lid for non-wetting crucible, Dia 56mm, fits Stock #s 46863 & 46458
46944	Platinum Lid for non-wetting crucible, Dia 59mm, fits Stock #s 46571 & 46804
46996	Platinum Lid for non-wetting crucible, Dia 62mm, fits Stock #s 46568 & 46429
46344	Platinum Lid for non-wetting crucible, Dia 65mm, fits Stock #s 46910 & 46850

Standard Non-Wetting Platinum 5% Gold Crucibles



High purity metals and materials
Fundamentals for research



[Discover our range >](#)








Platinum Labware




[Request a Quote >](#)

46376 Non-wetting platinum 5% gold crucible, Top Dia 26mm, Bot Dia 15.5mm, Ht 29mm, Base Thickness 0.25mm, Capacity 10mL

46483 Non-wetting platinum 5% gold crucible, Top Dia 29mm, Bot Dia 17.5mm, Ht 24mm, Base Thickness 0.25mm, Capacity 10mL

	46845	Non-wetting platinum 5% gold crucible, Top Dia 29mm, Bot Dia 17.5mm, Ht 33mm, Base Thickness 0.25mm, Capacity 15mL
	46900	Non-wetting platinum 5% gold crucible, Top Dia 32.5mm, Bot Dia 18mm, Ht 36mm, Base Thickness 0.25mm, Capacity 20mL
	46769	Non-wetting platinum 5% gold crucible, Top Dia 32mm, Bot Dia 19.5mm, Ht 27mm, Base Thickness 0.30mm, Capacity 15mL
	46859	Non-wetting platinum 5% gold crucible, Top Dia 35mm, Bot Dia 21mm, Ht 39mm, Base Thickness 0.30mm, Capacity 25mL
	46865	Non-wetting platinum 5% gold crucible, Top Dia 37.5mm, Bot Dia 22.5mm, Ht 41mm, Base Thickness 0.30mm, Capacity 30mL
	46822	Non-wetting platinum 5% gold crucible, Top Dia 39mm, Bot Dia 23.5mm, Ht 43mm, Base Thickness 0.30mm, Capacity 35mL
	46435	Non-wetting platinum 5% gold crucible, Top Dia 41.5mm, Bot Dia 25mm, Ht 46mm, Base Thickness 0.30mm, Capacity 40mL
	46542	Non-wetting platinum 5% gold crucible, Top Dia 41mm, Bot Dia 25mm, Ht 34mm, Base Thickness 0.30mm, Capacity 25mL
	46616	Non-wetting platinum 5% gold crucible, Top Dia 43.5mm, Bot Dia 26.5mm, Ht 36mm, Base Thickness 0.30mm, Capacity 35mL
	46701	Non-wetting platinum 5% gold crucible, Top Dia 44mm, Bot Dia 26.5mm, Ht 49mm, Base Thickness 0.38mm, Capacity 50mL
	46400	Non-wetting platinum 5% gold crucible, Top Dia 49mm, Bot Dia 30mm, Ht 41mm, Base Thickness 0.36mm, Capacity 50mL
	46901	Non-wetting platinum 5% gold crucible, Top Dia 50mm, Bot Dia 30mm, Ht 55mm, Base Thickness 0.41mm, Capacity 75mL
	46840	Non-wetting platinum 5% gold crucible, Top Dia 55mm, Bot Dia 33.5mm, Ht 46mm, Base Thickness 0.40mm, Capacity 70mL
	46863	Non-wetting platinum 5% gold crucible, Top Dia 56mm, Bot Dia 33mm, Ht 62mm, Base Thickness 0.48mm, Capacity 100mL
	46571	Non-wetting platinum 5% gold crucible, Top Dia 58mm, Bot Dia 35mm, Ht 64mm, Base Thickness 0.48mm, Capacity 125mL
	46568	Non-wetting platinum 5% gold crucible, Top Dia 62mm, Bot Dia 37mm, Ht 68mm, Base Thickness 0.48mm, Capacity 150mL
	46910	Non-wetting platinum 5% gold crucible, Top Dia 65mm, Bot Dia 39mm, Ht 71mm, Base Thickness 0.48mm, Capacity 175mL
	46438	Non-wetting platinum 5% gold crucible, Top Dia 68mm, Bot Dia 41mm, Ht 74mm, Base Thickness 0.48mm, Capacity 200mL

Platinum Crucible for Perl X 2 Machine with 5mm wide Flange




46663 Platinum crucible for Perl X 2 machine with 5mm wide flange, Top Dia 50mm, Bot Dia 50mm, Ht 40mm, Base Thickness 0.80mm



46496 Platinum crucible for Perl X 2 machine with 5mm wide flange, Top Dia 50mm, Bot Dia 50mm, Ht 40mm, Base Thickness 0.90mm


Platinum Crucible with Enlargement



46927


Platinum crucible with enlargement, Top Dia 34.5mm, Bot Dia 20mm, Ht 48mm, Base Thickness 0.40mm










Platinum Perl X Type Casting Moulds



46409	Platinum Leco type FX-100 & FX-200 Casting Mould, Top OD 44mm, Bot Dia 35mm, BT 0.76mm, Flange Shape: round, Ht 8mm
46955	Platinum Perl X type Casting Mould, Top OD 55mm, Top ID 32mm, Bot Dia 30.5mm, BT 0.5mm, Flange Shape: round, Ht 4mm
46978	Platinum Perl X type Casting Mould, Top OD 55mm, Top ID 32mm, Bot Dia 30.5mm, BT 0.7mm, Flange Shape: round, Ht 4mm
46404	Platinum Perl X type Casting Mould, Top OD 65mm, Top ID 31.5mm, Bot Dia 35.5mm, BT 1mm, Flange Shape: round, Ht 10mm
46373	Platinum Perl X type Casting Mould, Top OD 65mm, Top ID 40.5mm, Bot Dia 39.5mm, BT 1mm, Flange Shape: round, Ht 10mm
46550	Platinum Standard Casting Mould, Top OD 41mm, Top ID 31mm, Bot Dia 29mm, BT 0.8mm, Flange Shape: square, Ht 3mm
46528	Platinum Standard Casting Mould, Top OD 44mm, Top ID 36.5mm, Bot Dia 33.5mm, BT 0.8mm, Flange Shape: round, Ht 3mm
46491	Platinum Standard Casting Mould, Top OD 44mm, Top ID 36.5mm, Bot Dia 33.5mm, BT 1mm, Flange Shape: round, Ht 3mm
46757	Platinum Standard Casting Mould, Top OD 51mm, Top ID 36.5mm, Bot Dia 33.5mm, BT 0.8mm, Flange Shape: square, Ht 3mm
46414	Platinum Standard Casting Mould, Top OD 51mm, Top ID 36.5mm, Bot Dia 33.5mm, BT 1mm, Flange Shape: square, Ht 3mm
46676	Platinum Standard Casting Mould, Top OD 51mm, Top ID 41mm, Bot Dia 39mm, BT 0.8mm, Flange Shape: square, Ht 3mm
46929	Platinum Standard Casting Mould, Top OD 51mm, Top ID 41mm, Bot Dia 39mm, BT 1mm, Flange Shape: square, Ht 3mm

Platinum Combustion Boats



	46425	Platinum Combustion Boat, Top 12mm x 6mm, Bottom 11mm x 4mm, Height 4mm, Base Thickness 0.15mm, Capacity \approx 0.20mL
	46684	Platinum Combustion Boat, Top length 100mm, width 12mm, Ht 10mm, Base Thickness 0.20mm, Capacity \approx 11mL
	46625	Platinum Combustion Boat, Top length 10mm, width 4mm, Ht 3mm, Base Thickness 0.15mm, Capacity \approx 0.10mL
	46405	Platinum Combustion Boat, Top length 120mm, width 14mm, Ht 10mm, Base Thickness 0.20mm, Capacity \approx 15mL
	46862	Platinum Combustion Boat, Top length 140mm, width 12mm, Ht 12mm, Base Thickness 0.20mm, Capacity \approx 17.5mL
	46653	Platinum Combustion Boat, Top length 15mm, width 6mm, Ht 4mm, Base Thickness 0.15mm, Capacity \approx 0.15mL
	46464	Platinum Combustion Boat, Top length 160mm, width 12mm, Ht 14mm, Base Thickness 0.20mm, Capacity \approx 24mL
	46598	Platinum combustion boat, Top length 180mm, width 14mm, Ht 16mm, Base Thickness 0.20mm, Capacity \approx 36mL
	46473	Platinum Combustion Boat, Top length 200mm, width 16mm, Ht 18mm, Base Thickness 0.20mm, Capacity \approx 52.5mL

	46361	Platinum Combustion Boat, Top length 20mm, width 7mm, Ht 4mm, Base Thickness 0.15mm, Capacity ≈ 0.50mL
	46541	Platinum Combustion Boat, Top length 25mm, width 5mm, Ht 4mm, Base Thickness 0.15mm, Capacity ≈ 0.45mL
	46599	Platinum Combustion Boat, Top length 30mm, width 5mm, Ht 4mm, Base Thickness 0.15mm, Capacity ≈ 0.50mL
	46338	Platinum Combustion Boat, Top length 50mm, width 12mm, Ht 10mm, Base Thickness 0.15mm, Capacity ≈ 5.40mL
	46899	Platinum Combustion Boat, Top length 60mm, width 12mm, Ht 10mm, Base Thickness 0.15mm, Capacity ≈ 6.50mL
	46767	Platinum Combustion Boat, Top length 60mm, width 9mm, Ht 9mm, Capacity 3mL
	46449	Platinum Combustion Boat, Top length 70mm, width 9mm, Ht 8mm, Base Thickness 0.20mm, Capacity ≈ 4.5mL
	46508	Platinum Combustion Boat, Top length 80mm, width 10mm, Ht 8mm, Base Thickness 0.20mm, Capacity ≈ 6.0mL
	46573	Platinum Combustion Boat with handle, length 110mm, width 13mm, Ht 10mm, Capacity 12mL
	46818	Platinum Combustion Boat with handle, length 48mm, width 13mm, Ht 10mm, Capacity 4mL
	46713	Platinum Combustion Boat with handle, length 60mm, width 13mm, Ht 10mm, Capacity 6mL
	46891	Platinum Combustion Boat with handle, length 60mm, width 9mm, Ht 9mm, Capacity 3mL
	46648	Platinum Combustion Boat with handle, length 85mm, width 13mm, Ht 10mm, Capacity 9mL
	46363	Platinum Combustion Boat without handle, length 100mm, width 16mm, Ht 13mm, Capacity 12mL
	46524	Platinum Combustion Boat without handle, length 60mm, width 9mm, Ht 9mm, Capacity 3mL

Platinum Crucible for Schoeps Automatic Fusion Machine with Supporting Wire Cage and pouring lip



46546 Platinum crucible for Schoeps AF machine w/sup wire cage & pour lip, Top Dia 35mm, Bot Dia 21mm, Ht 39mm, BT 0.32mm



46463 Platinum crucible for Schoeps AF machine w/sup wire cage & pour lip, Top Dia 37.5mm, Bot Dia 22.5mm, Ht 40mm, BT 0.32mm



46357 Platinum crucible for Schoeps AF machine w/sup wire cage & pour lip, Top Dia 39mm, Bot Dia 23.5mm, Ht 43mm, BT 0.32mm















46697 Platinum crucible for Schoeps AF machine w/sup wire cage & pour lip, Top Dia 44mm, Bot Dia 26.8mm, Ht 49mm, BT 0.32mm


Platinum Electrodes



In constructing general service electrodes, platinum is the only material that offers sufficiently high electrical conductivity, mechanical strength and resistance to electrolyte and chemical attack. For most analysis, the best cathode is an ordinary gauze cylinder. This material offers greater surface area than equal weight of plain or perforated foil, and allows the electrolyte to circulate more freely.


	46707	Platinum Anode mesh cylinder, cyl. Dia 12mm, Ht 50mm, overall Ht 130mm
	46788	Platinum Anode mesh cylinder with central stem, cyl. Dia 25mm, Ht 45mm, overall Ht 170mm, ≈S.A. 70cm ²
	46921	Platinum Anode rotating paddle shaped, Ht 140mm, ≈S.A. 18cm ²
	46386	Platinum Anode wire helix, helix Dia 10mm, Ht 30mm, overall Ht 140mm, ≈S.A. 0.35cm ²
	46556	Platinum Anode wire helix, helix Dia 32mm, Ht 32mm, overall Ht 140mm, ≈S.A. 0.75cm ²
	46604	Platinum Cathode mesh cylinder, cyl. Dia 30mm, Ht 50mm, overall Ht 130mm
	46820	Platinum Cathode mesh cylinder, cyl. Dia 40mm, Ht 40mm, overall Ht 140mm, ≈S.A. 100cm ²
	46619	Platinum Fischer Anode mesh cylinder, cyl. Dia 32mm, Ht 32mm Ht 32mm, overall Ht 140mm, ≈S.A. 70cm ²
	46940	Platinum Fischer Cathode mesh cylinder with rings on stem, cyl. Dia 45mm, Ht 45mm, overall Ht 125mm, ≈S.A. 125cm ²
	46569	Platinum Fischer Cathode mesh cylinder w/out rings on stem, cyl. Dia 45mm, Ht 45mm, overall Ht 150mm, ≈S.A. 125cm ²
	46545	Platinum Shoniger Basket, Dia 7mm, Ht 20mm, overall Ht 70mm
	46396	Platinum Winkler Cathode divides mesh cylinder, cyl. Dia 21mm, Ht 45mm, overall Ht 150mm, ≈S.A. 120cm ²

Platinum Gauzes




40931	Gold gauze, 100 mesh woven from 0.064mm (0.0025in) dia wire, 99.99% (metals basis)
40930	Gold gauze, 52 mesh woven from 0.102mm (0.004in) dia wire, 99.99% (metals basis)
41814	Platinum gauze, 45 mesh woven from 0.198mm (0.0078in) dia wire, 99.9% (metals basis)
10283	Platinum gauze, 52 mesh woven from 0.1mm (0.004in) dia wire, 99.9% (metals basis)
46966	Platinum gauze, Unimesh N7433, expanded metal mesh, 0.34mm (0.013in) thick
40935	Silver gauze, 50 mesh woven from 0.0764mm (0.003in) dia wire

Platinum Inoculation Loops




20965	Pt 15% Ir Inoculating Loop, Loop ID 2mm, Wire Dia 0.41mm, Shank length 45mm, Delivery Cap. 0.001mL, ID +0.01/-0.06mm
21227	Pt 15% Ir Inoculating Loop, Loop ID 3mm, Wire Dia 0.51mm, Shank length 45mm, Delivery Cap. 0.003mL, ID +0.01/-0.06mm
21233	Pt 3.5% Rh Inoculating Loop, Loop ID 1.45mm, Wire Dia 0.41mm, Shank length 25mm, Delivery Cap. 0.001mL, ID +0.01/-0.06mm
98123	Pt 3.5% Rh Inoculating Loop, Loop ID 1.45mm, Wire Dia 0.41mm, Shank length 75mm, Delivery Cap. 0.001mL, ID +0.01/-0.06mm
21231	Pt 3.5% Rh Inoculating Loop, loop ID 4mm, Wire Dia 0.91mm, Shank length 25mm, Delivery Cap. 0.011mL, ID +0.01/-0.06mm
21236	Pt 3.5% Rh Inoculating Loop, Loop ID 4mm, Wire Dia 0.91mm, Shank length 75mm, Delivery Cap. 0.011mL, ID +0.01/-0.06mm
21235	Pt 5% Ru Inoculating Loop, Loop ID 3mm, Wire Dia 0.41mm, Shank length 45mm, Delivery Cap. 0.003mL, ID +0.01/-0.06mm

Platinum Square Edge Dish



46778 Platinum Square Edge Dish, Bot Dia 40.9mm, Capacity 3.5mL

Platinum Crucible for Herzog Automatic Fusion Machine



46949 Platinum crucible for Herzog automatic fusion machine, Top Dia 50mm, Bot Dia 38.5mm, Ht 30mm, Base Thickness 1.0mm


Platinum Crucible for Schoeps Semi-Automatic Fusion Machine with Supporting Pins



46874 Platinum crucible for Schoeps semi-auto fusion machine w/sup pins, Top Dia 40mm, Bot Dia 30mm, Ht 30mm, Bt 0.45mm

46947 Platinum crucible for Schoeps semi-auto fusion machine w/sup pins, Top Dia 40mm, Bot Dia 30mm, Ht 30mm, Bt 0.60mm

Platinum Evaporating Dishes



Alfa Aesar evaporating dishes can be supplied with reinforced tops or bottoms as well as covers. Most dishes are constructed of platinum but other alloys are available. Alfa Aesar can also assist with custom fabricated dishes to your own designs and specifications.


Platinum Classen Dish



46826

Platinum Classen dish, Bot Dia 90mm, Ht 45mm, Base Thickness 0.15mm, Capacity 250mL

Platinum Round Bottom Dishes with Handles



46564


Platinum Round Bottom Dish with Handles, Top Dia 15mm, Ht 8mm











46957


Platinum Round Bottom Dish with Handles, Top Dia 22mm, Ht 11mm






Platinum Covers for Standard Dishes




	46915	Platinum Lid for Standard Dish, Dia 35mm, fits stock #s 46987 & 46893
	46431	Platinum Lid for Standard Dish, Dia 38mm, fits stock #s 46878 & 46642
	46844	Platinum Lid for Standard Dish, Dia 42mm, fits stock #s 46729 & 46821
	46456	Platinum Lid for Standard Dish, Dia 50mm, fits stock #s 46932 & 46680
	46582	Platinum Lid for Standard Dish, Dia 57.5mm, fits stock #s 46775 & 46424
	46621	Platinum Lid for Standard Dish, Dia 62.5mm, fits stock #s 46994 & 46842
	46981	Platinum Lid for Standard Dish, Dia 70mm, fits stock #s 46611 & 46742
	46937	Platinum Lid for Standard Dish, Dia 75mm, fits stock #s 46711 & 46919

Platinum Standard Dishes



	46352	Platinum Standard Dish, Top Dia 100mm, Ht 41mm, Base Thickness 0.33mm, Capacity 250mL
	46576	Platinum Standard Dish, Top Dia 125mm, Ht 50mm, Base Thickness 0.32mm, Capacity 500mL
	46397	Platinum Standard Dish, Top Dia 143mm, Ht 58mm, Base Thickness 0.33mm, Capacity 700mL
	46987	Platinum Standard Dish, Top Dia 35mm, Ht 17mm, Base Thickness 0.17mm, Capacity 10mL
	46878	Platinum Standard Dish, Top Dia 38mm, Ht 19mm, Base Thickness 0.17mm, Capacity 15mL
	46729	Platinum Standard Dish, Top Dia 42mm, Ht 20mm, Base Thickness 0.17mm, Capacity 20mL
	46932	Platinum Standard Dish, Top Dia 50mm, Ht 22mm, Base Thickness 0.17mm, Capacity 30mL
	46775	Platinum Standard Dish, Top Dia 57.5mm, Ht 26mm, Base Thickness 0.17mm, Capacity 50mL
	46994	Platinum Standard Dish, Top Dia 62.5mm, Ht 27mm, Base Thickness 0.17mm, Capacity 60mL
	46611	Platinum Standard Dish, Top Dia 70mm, Ht 28mm, Base Thickness 0.17mm, Capacity 80mL
	46711	Platinum Standard Dish, Top Dia 75mm, Ht 29mm, Base Thickness 0.23mm, Capacity 100mL
	46399	Platinum Standard Dish, Top Dia 82.5mm, Ht 34mm, Base Thickness 0.25mm, Capacity 150mL
	46782	Platinum Standard Dish, Top Dia 90mm, Ht 38mm, Base Thickness 0.25mm, Capacity 175mL
	46905	Platinum Standard Dish, Top Dia 90mm, Ht 40mm, Base Thickness 0.28mm, Capacity 200mL

Platinum Flat Bottom Plain Dishes



46889 Platinum Flat Bottom Plain Dish, Top Dia 10mm, Ht 2mm




46590 Platinum Flat Bottom Plain Dish, Top Dia 15mm, Ht 2mm

Platinum Standard Dishes with Reinforced Rims

	46928	Platinum Standard Dish with reinforced rim, Top Dia 100mm, Ht 41mm, Base Thickness 0.33mm, Capacity 250mL
	46614	Platinum Standard Dish with reinforced rim, Top Dia 125mm, Ht 50mm, Base Thickness 0.33mm, Capacity 500mL
	46689	Platinum Standard Dish with reinforced rim, Top Dia 143mm, Ht 58mm, Base Thickness 0.33mm, Capacity 700mL
	46893	Platinum Standard Dish with reinforced rim, Top Dia 35mm, Ht 17mm, Base Thickness 0.17mm, Capacity 10mL
	46642	Platinum Standard Dish with reinforced rim, Top Dia 38mm, Ht 19mm, Base Thickness 0.17mm, Capacity 15mL
	46821	Platinum Standard Dish with reinforced rim, Top Dia 42mm, Ht 20mm, Base Thickness 0.17mm, Capacity 20mL
	46680	Platinum Standard Dish with reinforced rim, Top Dia 50mm, Ht 22mm, Base Thickness 0.17mm, Capacity 30mL
	46424	Platinum Standard Dish with reinforced rim, Top Dia 57.5mm, Ht 26mm, Base Thickness 0.17mm, Capacity 50mL
	46842	Platinum Standard Dish with reinforced rim, Top Dia 62.5mm, Ht 27mm, Base Thickness 0.17mm, Capacity 60mL
	46742	Platinum Standard Dish with reinforced rim, Top Dia 70mm, Ht 28mm, Base Thickness 0.17mm, Capacity 80mL
	46919	Platinum Standard Dish with reinforced rim, Top Dia 75mm, Ht 29mm, Base Thickness 0.23mm, Capacity 100mL
	46457	Platinum Standard Dish with reinforced rim, Top Dia 82.5mm, Ht 31mm, Base Thickness 0.23mm, Capacity 125mL
	46953	Platinum Standard Dish with reinforced rim, Top Dia 82.5mm, Ht 34mm, Base Thickness 0.25mm, Capacity 150mL
	46507	Platinum Standard Dish with reinforced rim, Top Dia 90mm, Ht 38mm, Base Thickness 0.25mm, Capacity 175mL
	46807	Platinum Standard Dish with reinforced rim, Top Dia 90mm, Ht 40mm, Base Thickness 0.28mm, Capacity 200mL

Platinum Round Bottom Dishes



46461 Platinum Round Bottom Dish, Top OD 100mm, Ht 50mm, Capacity 250mL



46387 Platinum Round Bottom Dish, Top OD 110mm, Ht 55mm, Capacity 325mL



46674 Platinum Round Bottom Dish, Top OD 120mm, Ht 60mm, Capacity 400mL



46753 Platinum Round Bottom Dish, Top OD 130mm, Ht 65mm, Capacity 500mL



46898 Platinum Round Bottom Dish, Top OD 40mm, Ht 20mm, Capacity 20mL



46761 Platinum Round Bottom Dish, Top OD 50mm, Ht 25mm, Capacity 35mL



46839 Platinum Round Bottom Dish, Top OD 60mm, Ht 30mm, Capacity 50mL



46570 Platinum Round Bottom Dish, Top OD 70mm, Ht 35mm, Capacity 90mL




46385 Platinum Round Bottom Dish, Top OD 80mm, Ht 40mm, Capacity 150mL



46471 Platinum Round Bottom Dish, Top OD 90mm, Ht 45mm, Capacity 200mL

Platinum Wine Dish



46512 Platinum wine dish, Bot Dia 85mm, Ht 20mm, Base Thickness 0.15mm

Platinum Tongs & Tweezers




46362 Bowed Stainless Steel Tong with Platinum Tips, Length 200mm




46638 Stainless Steel Tweezer with Platinum Tips, Length 130mm

Platinum Adjustable Triangles




	46802	Platinum adjustable triangle, Width 35mm, Length 60mm, Base Thickness 1.2mm
	46858	Platinum adjustable triangle, Width 40mm, Length 88mm, Base Thickness 1.2mm
	46668	Platinum adjustable triangle, Width 50mm, Length 95mm, Base Thickness 1.2mm
	46445	Platinum adjustable triangle, Width 60mm, Length 115mm, Base Thickness 1.5mm
	46967	Platinum adjustable triangle, Width 70mm, Length 120mm, Base Thickness 1.5mm
	46662	Platinum adjustable triangle, Width 80mm, Length 148mm, Base Thickness 1.5mm

Platinum Fixed Triangles



	46358	Platinum fixed triangle, Width 30mm, Length 20mm, Base Thickness 1.2mm
	46440	Platinum fixed triangle, Width 35mm, Length 20mm, Base Thickness 1.2mm
	46394	Platinum fixed triangle, Width 40mm, Length 20mm, Base Thickness 1.2mm
	46933	Platinum fixed triangle, Width 50mm, Length 25mm, Base Thickness 1.2mm
	46917	Platinum fixed triangle, Width 60mm, Length 25mm, Base Thickness 1.2mm
	46759	Platinum fixed triangle, Width 70mm, Length 25mm, Base Thickness 1.2mm
	46876	Platinum fixed triangle, Width 80mm, Length 25mm, Base Thickness 1.2mm
	46733	Platinum fixed triangle, Width 90mm, Length 25mm, Base Thickness 1.2mm

Platinum Shovels



	46852	Platinum shovel, Width 10mm, Length 50mm, Base Thickness 1.5mm
	46885	Platinum shovel, Width 15mm, Length 60mm, Base Thickness 1.5mm
	46704	Platinum shovel, Width 20mm, Length 80mm, Base Thickness 1.5mm
	46817	Platinum shovel, Width 25mm, Length 100mm, Base Thickness 2mm
	46600	Platinum shovel, Width 25mm, Length 50mm, Base Thickness 2mm
	46679	Platinum shovel, Width 5mm, Length 50mm, Base Thickness 1mm

Platinum Stirring Spatulas

	46777	Platinum stirring spatula double ended, Width 10mm, Length 100mm, Base Thickness 2mm
	46706	Platinum stirring spatula double ended, Width 10mm, Length 150mm, Base Thickness 2mm
	46608	Platinum stirring spatula double ended, Width 5mm, Length 50mm, Base Thickness 1mm
	46340	Platinum stirring spatula double ended, Width 7mm, Length 75mm, Base Thickness 2mm
	46617	Platinum stirring spatula single ended, Width 10mm, Length 100mm, Base Thickness 2mm
	46560	Platinum stirring spatula single ended, Width 10mm, Length 150mm, Base Thickness 2mm
	46765	Platinum stirring spatula single ended, Width 5mm, Length 50mm, Base Thickness 1mm
	46612	Platinum stirring spatula single ended, Width 7mm, Length 75mm, Base Thickness 2mm

Алматы (7273)495-231
 Ангарск (3955)60-70-56
 Архангельск (8182)63-90-72
 Астрахань (8512)99-46-04
 Барнаул (3852)73-04-60
 Белгород (4722)40-23-64
 Благовещенск (4162)22-76-07
 Брянск (4832)59-03-52
 Владивосток (423)249-28-31
 Владикавказ (8672)28-90-48
 Владимир (4922)49-43-18
 Волгоград (844)278-03-48
 Вологда (8172)26-41-59
 Воронеж (473)204-51-73
 Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
 Ижевск (3412)26-03-58
 Иркутск (395)279-98-46
 Казань (843)206-01-48
 Калининград (4012)72-03-81
 Калуга (4842)92-23-67
 Кемерово (3842)65-04-62
 Киров (8332)68-02-04
 Коломна (4966)23-41-49
 Кострома (4942)77-07-48
 Краснодар (861)203-40-90
 Красноярск (391)204-63-61
 Курск (4712)77-13-04
 Курган (3522)50-90-47
 Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
 Москва (495)268-04-70
 Мурманск (8152)59-64-93
 Набережные Челны (8552)20-53-41
 Нижний Новгород (831)429-08-12
 Новокузнецк (3843)20-46-81
 Ноябрьск (3496)41-32-12
 Новосибирск (383)227-86-73
 Омск (3812)21-46-40
 Орел (4862)44-53-42
 Оренбург (3532)37-68-04
 Пенза (8412)22-31-16
 Петрозаводск (8142)55-98-37
 Псков (8112)59-10-37
 Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
 Рязань (4912)46-61-64
 Самара (846)206-03-16
 Санкт-Петербург (812)309-46-40
 Саратов (845)249-38-78
 Севастополь (8692)22-31-93
 Саранск (8342)22-96-24
 Симферополь (3652)67-13-56
 Смоленск (4812)29-41-54
 Сочи (862)225-72-31
 Ставрополь (8652)20-65-13
 Сургут (3462)77-98-35
 Сыктывкар (8212)25-95-17
 Тамбов (4752)50-40-97
 Тверь (4822)63-31-35

Тольятти (8482)63-91-07
 Томск (3822)98-41-53
 Тула (4872)33-79-87
 Тюмень (3452)66-21-18
 Ульяновск (8422)24-23-59
 Улан-Уде (3012)59-97-51
 Уфа (347)229-48-12
 Хабаровск (4212)92-98-04
 Чебоксары (8352)28-53-07
 Челябинск (351)202-03-61
 Череповец (8202)49-02-64
 Чита (3022)38-34-83
 Якутск (4112)23-90-97
 Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(7172)727-132

Киргизия +996(312)96-26-47

<https://aesar.nt-rt.ru/> || arj@nt-rt.ru